


BIOVIA PIPELINE PILOT SEQUENCE ANALYSIS COMPONENT COLLECTION

DATASHEET

The BIOVIA Pipeline Pilot Sequence Analysis offers essential bioinformatics tools and algorithms for creating practical sequence analysis workflows that complement any laboratory research. With over 100 different component functions, you can analyze and annotate DNA and protein sequences using a variety of industry accepted methods.

WITH THE BIOVIA PIPELINE PILOT SEQUENCE ANALYSIS, YOU CAN:

- Automate your routine analysis of new sequence data
- Compare genes and genomes across species
- Prototype new algorithmic processes
- Perform chemo-genomic analysis with the Chemistry Collection
- Integrate your own algorithms into any protocol



BIOVIA Pipeline Pilot Sequence Analysis can be combined into logical networks that automate the analysis and assembly of sequence families.

SEQUENCE SIMILARITY SEARCHING

For uncomplicated sequence similarity searching, use the Smith-Waterman component or any of our standard BLAST components, including BLASTn, BLASTp, BLASTx, tBLASTn, and tBLASTx. Display the results of these searches with our Similarity Search Viewer, and extract or fetch to further analyze the individual hits. BIOVIA Pipeline Pilot also makes it easy to build custom BLAST databases on-the-fly, providing more targeted and intelligent sequence similarity searches.

MANIPULATION AND ANNOTATION

An assortment of sequence manipulation and annotation components are available. For DNA sequences, these components include functions such as primer identification, GC content, six-frame translations, reverse complement, and siRNA target site prediction. For protein sequences, you can determine back translation, secondary structure prediction, and isoelectric point.

ALIGNMENT AND PROFILE SEARCHING

To identify potential homologs among a collection of different organisms and quickly select variants among nucleotide sequence regions, you can use our alignment and profiling components. To allow for multiple sequence alignments of either DNA or protein sequences, ClustalW is included. For profiling tasks, HMM Build, Align, Search, and Pfam are available.

PATTERN MATCHING

A variety of standard tools are implemented so you can search for interesting patterns or motifs within a biological sequence.

These algorithms enable the identification of potential PROSITE regions, GC rich regions, proteolytic cleavage sites, restriction enzyme sites, signal peptide cleavage sites, open reading frames, or regular expression patterns.

RESULT VIEWERS

You can also visualize sequence information and related features with our sequence viewers. These viewers include a custom report (PDF or HTML), plain text view, or Artemis viewer (for nucleotide sequences). Also, multiple sequence alignments can be displayed in a custom report (PDF or HTML), plain text, or lalView.

THIRD-PARTY TOOL INTEGRATION

The BIOVIA Pipeline Pilot Sequence Analysis includes examples of integration with BioPerl, NCBI BLAST, GCG programs, EMBOSS tools, and BioJava. You can use these examples as templates to extend the available functionality to include other programs of interest. Your BIOVIA Pipeline Pilot integration options include Java, Perl, SOAP,VB Script, or writing simple command-line wrappers.

To learn more about BIOVIA Pipeline Pilot, go to accelrys.com/pipeline-pilot

Our **3D**EXPERIENCE Platform powers our brand applications, serving 12 industries, and provides a rich portfolio of industry solution experiences.

Dassault Systèmes, the **3D**EXPERIENCE Company, provides business and people with virtual universes to imagine sustainable innovations. Its world-leading solutions transform the way products are designed, produced, and supported. Dassault Systèmes' collaborative solutions foster social innovation, expanding possibilities for the virtual world to improve the real world. The group brings value to over 170,000 customers of all sizes in all industries in more than 140 countries. For more information, visit **www.3ds.com**.



3DEXPERIENCE

Dassault Systèmes Corporate Dassault Systèmes 175 Wyman Street Waltham, Massachusetts 02451-1223 BIOVIA Corporate Americas BIOVIA 5005 Wateridge Vista Drive, San Diego, CA 92121 BIOVIA Corporate Europe BIOVIA 334 Cambridge Science Park, Cambridge CB4 0WN England